




# Welcome to Pretaa!

Let's take a look at how we work for you.



Once you're signed in, you'll begin with Pretaa's Clinician's view on the **Events** tab.

This is where Significant Patient Events are shown such as Daily Reports, Anomalies, and Assessment results Filter as necessary for a specific alert or status. Drill down into the details of the specific event right from here.

**Events** [Search] [Date Range: 08/18/2023 - 08/18/2023]

- ALERT **Breanr** [ ] ALERT 08/05/23 | ...
- [ ] COMPLETED ASSESSMENT
- ALERT **Jaylon** [ ] REPORT 08/05/23 | ...
- [ ] CONTACTED HELPLINE
- ALERT **Jaiden** [ ] FENCE 08/05/23 | ...
- [ ] IN CENSUS
- REPORT **Daily F** [ ] DISCHARGED 08/05/23 | ...

**Event Details**

ALERT **Jaiden Cremin's** data shows anomalous HRV readings. 08/06/23

**DAILY REPORT**  
Event date: 08/03/2023  
Date published: 08/04/2023  
Jaiden Cremin  
Patient Report: 8286167668740803

Biometric Scale

- Heart Rate (BPM)
- Steps (Daily Total)
- Heart Rate Variability (Milliseconds)
- Temp (°F)

**Assessment details** 08/05/2023

COMPLETED ASSESSMENT  
URICA: University of Rhode Island Change Assessment was submitted on 08/05/2023 by Orpha Cassin

Participant Name: Orpha Cassin  
Initiated by: Snehl Adhikari  
Date completed: 08/05/2023

Readiness Scoring: 9.00

Score	Result
<8	Precontemplation
8-11	Precontemplation/Contemplation Transition
11-14	Contemplation
>11	Contemplation/Preparation Transition
>14	Prep/Action

Please indicate the extent to which you tend to agree or disagree with each statement. In each case, make your choice in terms of how you feel right now, and what you know *is* in the past or would like to feel. There are CIVIC possible responses to each of the items in the questionnaire.

On the **Patient's** tab, you will find everyone you are responsible for and be able to filter and select only those needing review. From there, you can drill down into each patient and view their details, any notes you have added, their contacts, their timeline of all events, assessments, and any related geofences.

**Patients List Screenshot:**

- Search bar: Q Search
- Filter modal:
  - Clear All
  - Eyeball highlighted
  - Eyeball not highlighted
  - In Census
  - Discharged
  - Submit
- Table of patients:
 

Group	Name	Eye Icon
A	Aliza Jenkins	👁
	Anthony Alvarez	👁
B	Bessie Johnston	👁
	Billy Mills	👁
D		👁
		👁
		👁
		👁
		👁
		👁
		👁

**Patient Detail Screenshot (Jaylon Murray):**

- Location: Tranquil Detox Retreat
- Emergency Contact: Liam Murray
- Age: 33
- Payment Method: Insurance
- Diagnosis: N/A
- Active: YES
- Intake Date: 05/16/2023
- Bed Name: Bed 1
- Insurance company: A.T.C.
- Solitary Days: 94
- Discharge Date: 06/27/2023
- Room #: Pine
- Payment method Category: -

**Timeline Screenshot (Jaylon Murray's Timeline):**

- Event 1: FINCE - Jaylon Murray Has Entered Skid Row At 07:10 PM. (08/06/23)
- Event 2: REPORT - Jaylon Murray's Daily Report Is Now Available. (08/06/23)
- Event 3: ALERT - Jaylon Murray's Data Shows Anomalous Heart Readings. (08/05/23)
- Event 4: REPORT - Jaylon Murray's Daily Report Is Now Available. (08/05/23)
- Event 5: REPORT - Jaylon Murray's Daily Report Is Now Available. (08/04/23)
- Event 6: REPORT - Jaylon Murray's Daily Report Is Now Available. (08/03/23)
- Event 7: REPORT - Jaylon Murray's Daily Report Is Now Available. (08/02/23)

The **Assessments** tab is where you will find the clinically validated assessments. These include the BAM R and IOP, the PHQ-9 and 15, GAD-7, and URICA. These can be set up as a one time event or on a recurring basis, by clicking on the desired assessment, and clicking schedule new. There is also a place to schedule custom assessments that have been created. These can be sent to a single patient or any group of patients.

- Events
- Patients
- Assessments
- Reporting
- Settings
- Feedback
- Logout
- Contact Pretaa
- Version 1.1.0

## Assessment Templates

Clinical
Custom

Template Name	Topic	# of Assessments	
<b>PHQ-9 Assessment</b> Patient Depression Questionnaire (PHQ-9)	Depression	0	...
<b>PHQ-15 Assessment</b> The PHQ-15 is a brief, self-administered questionnaire used to assess somatic symptom severity and the potential presence of somatization and somatoform disorders.	Somatization	0	...
<b>GAD-7 Assessment</b> GAD-7 Generalized Anxiety Disorder	Anxiety	0	...
<b>BAM-R Assessment</b> Brief addiction monitor assessment (BAM) for continuous response (R)	Health, Alcohol, and Drug use	0	...
<b>BAM-IOP Assessment</b> Brief addiction monitor assessment (BAM) for intensive outpatient treatment (IOP)	Health, Alcohol, and Drug use	0	...
<b>URICA: University of Rhode Island Change Assessment</b> URICA: University of Rhode Island Change Assessment	Readiness	0	...

- Events
- Patients
- Assessments
- Reporting
- Settings
- Feedback
- Logout
- Contact Pretaa
- Version 1.1.0

< BACK
Template Name  
PHQ-9 Assessment

## Select Patients

Akeem Hickie x
Alda Williamson x

<input type="checkbox"/>	First Name	Last Name	Email	Care Teams	Patient Status	Gender
<input checked="" type="checkbox"/>	Akeem	Hickie	julio_heidenreich@yah...	N/A	Discharged	Male
<input checked="" type="checkbox"/>	Alda	Williamson	sallie53@gmail.com	N/A	Discharged	Female
<input type="checkbox"/>	Alessandra	Walter	karine.dietrich@yaho...	N/A	Discharged	Female
<input type="checkbox"/>	Alexane	Bradtkie	hazel_ohara5@hotmail...	N/A	Discharged	Female
<input type="checkbox"/>	Ali	Gutkowski	norwood_nitzsche84@...	N/A	Discharged	Male
<input type="checkbox"/>	Alison	Lind	dock44@hotmail.com	N/A	Discharged	Female
<input type="checkbox"/>	Aliza	Jenkins	frederick34@yahoo.com	N/A	In	Female
<input type="checkbox"/>	Annie	Swift	edison.adams@yahoo...	N/A	Discharged	Male

Rows: 123

**Select assessment date**

Now OR 
 Signature Required

Send
Cancel


The **Reporting** tab includes the **Events** tab where you can review anomalies reported, assessments needing attention, suicidal ideations, helpline contacts, and geofences compromised. Filters allow for research by time or by patient(s). Individual patient information can also be accessed from this page by clicking on the patient's name.

- Events
- Patients
- Assessments
- Reporting**
  - Events
  - Assessment Stats
- Settings
- Feedback
- Logout


## Reporting/Events

Filter By: 3 Months | All patients selected

**322** Anomalies Reported | **263** Assessment Needing Attention | **43** Suicidal Ideation | **2** Help Line Contacted | **4** Geofences Compromised



- Events
- Patients
- Assessments
- Reporting**
  - Events
  - Assessment Stats
- Settings
- Feedback
- Logout



Patients

Name	Total	Biometric Scale
Dejuan Moen	1	🔴🟡🟢🟣
Nicolas O'Connell	1	🔴🟡🟢🟣
Bonnie Predovic	11	🔴🟡🟢🟣
Eddi Doyle	2	🔴🟡🟢🟣
Aliza Jenkins	0	🔴🟡🟢🟣
Donato Franekki	5	🔴🟡🟢🟣
Hulda Schmeier	6	🔴🟡🟢🟣

Row: 122

- Events
- Patients
- Reporting**
  - Events
  - Assessment Stats
- Settings
- Feedback
- Logout

### Reporting/Events - Jaylon Murray

Filter By: Month To Date

**2** Anomalies Reported | **0** Assessment Needing Attention | **0** Suicidal Ideation | **0** Help Line Contacted | **1** Geofences Compromised



**Assessment Stats** on the **Reporting** tab will provide statistics from the results of clinically validated assessments that the patients have self-administered. You can view these filtered by time, by patient, or by assessment type. Performance from admission to discharge by type of assessment can also be reviewed here. Results by severity by assessment type can be reviewed in the details by clicking on the desired assessment.

**pretaa**

Events

Patients

Assessments

**Reporting**

Events

[Assessment Stats](#)

Settings

Feedback

Logout

Contact Pretaa

### Reporting: Assessments

Filter By: 3 Months | All | All patients selected

Patients Overview	<b>309</b> PHQ-9	<b>167</b> BAM-R	<b>319</b> GAD-7	<b>303</b> URICA
-------------------	---------------------	---------------------	---------------------	---------------------

#### Number of patient Assessments by type

Assessment	PHQ-9	GAD-7	URICA	BAM-R	INCOMPLETE
First Assessment	20	20	20	20	5
Second Assessment	15	20	20	20	25
Third Assessment	15	20	20	20	25
Fourth Assessment	15	20	20	20	25
Fifth Assessment	15	20	20	20	25
Sixth Assessment	10	15	15	15	35

**pretaa**

Events

Patients

**Reporting**

Events

[Assessment Stats](#)

Settings

Feedback

Logout

Contact Pretaa

Version 1.1.0

#### Performance from admission to discharge

PHQ-9	↓ 55%	Average Score
GAD-7	↓ 49%	Average Score
URICA	↓ 23%	Readiness Score
BAM-R	↓ 99%	Use Factors

**pretaa**

Events

Patients

Assessments

**Reporting**

Events

[Assessment Stats](#)

Settings

Feedback

Logout

Contact Pretaa

Version 1.1.0

### Reporting: Assessments

Filter By: 3 Months | All | All patients selected

Patients Overview	<b>309</b> PHQ-9	<b>167</b> BAM-R	<b>319</b> GAD-7	<b>303</b> URICA
-------------------	---------------------	---------------------	---------------------	---------------------

#### PHQ-9 results by severity

Assessment	Severe	Moderately severe	Moderate	Mild	None	Average Score
First PHQ-9	18	2	2	2	2	20
Second PHQ-9	15	2	2	2	2	15
Third PHQ-9	10	2	2	2	2	10
Fourth PHQ-9	10	2	2	2	2	10
Fifth PHQ-9	10	2	2	2	2	10
Sixth PHQ-9	10	2	2	2	2	10

5

Welcome to Pretaa! Let's take a look at how we work for you.

Once you're signed in, you'll begin with Pretaa's Clinician's view on the Events tab. This is where Significant Patient Events are shown such as Daily Reports, Anomalies, and Assessment results Filter as necessary for a specific alert or status. Drill down into the events or assessment results right from here.

The collage displays three main screenshots from the Pretaa interface:

- Events Tab:** Shows a list of events with a search bar and a date range filter (08/18/2023 - 08/18/2023). A filter dropdown menu is open, listing event types: ALERT, COMPLETED ASSESSMENT, REPORT, CONTACTED HELPLINE, FENCE, IN CENSUS, and DISCHARGED. An 'Apply' button is at the bottom of the filter.
- Event Details:** Shows details for an alert titled "Jaidden Cremin's data shows anomalous HRV readings." It includes a "DAILY REPORT" section with a "Biometric Scale" gauge and three horizontal progress bars for:
  - Heart Rate (BPM): 71
  - Steps (Daily total): 7.7K
  - Heart Rate Variability (Milliseconds): 21
- Assessment details:** Shows details for an assessment titled "University of Rhode Island Change Assessment was submitted on 08/05/2023 by Orpha Cassin". It includes a "Patient Name" field, a "Result: Contemplation" status, a score of 9.00, a line graph showing scores over time, and a table of score ranges and corresponding transition types.
 

Score	Result: Contemplation
<8	Precontemplation
8-11	Precontemplation/Contemplation Transition
8-11	Contemplation
-11	Contemplation/Preparation Transition
11-14	Prep/Action

On the Patient's tab, you will find everyone you are responsible for and be able to filter and view only those needing review. From there, you can drill down into each patient and view their details, any notes you have added, their contacts, their timeline of all events, assessments, and any related geofences.

The screenshot shows the 'Patients' tab in the Pretaa application. A search bar is at the top with a placeholder 'Q Search'. Below it is a button labeled 'In Census X'. A modal window titled 'Filter' is open, showing options to filter patients: 'Clear All', 'Eyeball highlighted', 'Eyeball not highlighted', 'In Census' (checked), and 'Discharged'. A yellow 'Submit' button is at the bottom of the modal. The background shows a list of patients grouped by letter: A (Aliza Jenkins, Anthony Alvarez), B (Bessie Johnston, Billy Mills), D (Dorcas Simonis), E (Ellsworth Hahn, Emmitt Ebert). Each patient name has an eyeball icon to its right. A 'Clear all' link is in the top right corner of the patient list area. A dark overlay at the top right says 'Press F11 to exit full screen'.

The screenshot shows the patient details page for 'Jaylon Murray'. The page is divided into several sections:
 

- Header:** 'Jaylon Murray' with a '< BACK' link and an 'Invite Supporter' button.
- Location:** Tranquill Detox Retreat. Status: Active YES.
- Emergency Contact:** Liam Murray. Intake Date: 05/16/2023.
- Age:** 33. Bed Name: Bed 1.
- Payment Method:** Insurance company.
- Insurance:** N/A.
- Diagnosis:** N/A. Gender: Male.
- Notes:** Section for adding notes.
- Contacts:** Section for adding contacts.
- Timeline (176):** Section for viewing the patient's timeline.
- Assessments:** Section for viewing assessments.
- Geofences:** Section for viewing geofences.
- Last Known Location:** Section for viewing the last known location.

 On the right side, there is a 'Jaylon Murray's Timeline' section with a search bar and a list of events:
 

- FENCE:** Jaylon Murray Has Entered Skid Row At 07:10 PM. (08/06/23)
- REPORT:** Jaylon Murray's Daily Report Is Now Available. (08/06/23)
- ALERT:** Jaylon Murray's Data Shows Anomalous Heart Readings. (08/05/23)
- REPORT:** Jaylon Murray's Daily Report Is Now Available. (08/05/23)
- REPORT:** Jaylon Murray's Daily Report Is Now Available. (08/04/23)
- REPORT:** Jaylon Murray's Daily Report Is Now Available. (08/03/23)
- REPORT:** Jaylon Murray's Daily Report Is Now Available. (08/03/23)
- REPORT:** Jaylon Murray's Daily Report Is Now Available. (08/02/23)

 The left sidebar contains navigation options: Events, Patients (selected), Reporting, Settings, Feedback, Logout, Contact Pretaa, and Version 1.1.0.

The Reporting tab includes the Events tab where you can review anomalies reported, assessments needing attention, suicidal ideations, helpline contacts, and geofences compromised. Filters allow for research by time or by patient(s). Individual patient information can also be accessed from this page by clicking on the patient's name.

**pretaa**

Events

Patients

**Reporting**

Events

Assessment Stats

Settings

Feedback

Logout

Contact Pretaa

Version 1.1.0

## Reporting/Events

Filter By: 1 Month All patients selected

**160** Anomalies Reported

**138** Assessment Needing Attention

**22** Suicidal Ideation

**0** Help Line Contacted

**1** Geofences Compromised



### Types Of Anomalies

**pretaa**

Events

Patients

**Reporting**

Events

Assessment Stats

Settings

Feedback

Logout

Contact Pretaa

Version 1.1.0

Types Of Anomalies

Press **F11** to exit full screen

**160** Anomalies Reported

- 32 HRV Anomaly
- 55 Heart Anomaly
- 35 Sleep Anomaly
- 20 SpO2 Anomaly
- 18 Temp Anomaly

Name	Total	Biometric Scale	HRV Anomaly	Heart Anon
Ben Cleigley	1	<span style="color: red;">●</span>	0	0
Jaiden Cremin	6	<span style="color: red;">●</span>	2	0
Brenda Johnson	11	<span style="color: red;">●</span>	2	3
Fanny Nitzsche	1	<span style="color: red;">●</span>	0	1
Wolke Romaguera	11	<span style="color: red;">●</span>	3	0
Edd Doyle	1	<span style="color: red;">●</span>	1	0
Anthony Alvarez	3	<span style="color: red;">●</span>	1	0

Rows: 6/8

**pretaa**

Events

Patients

**Reporting**

Events

Assessment Stats

Settings

Feedback

Logout

Contact Pretaa

Version 1.1.0

Reporting/Events - Jaylon Murray

Filter By: Month To Date

**2** Anomalies Reported

**0** Assessment Needing Attention

**0** Suicidal Ideation

**0** Help Line Contacted

**1** Geofences Compromised

Date Range	HRV Anomaly	Heart Anomaly	Sleep Anomaly	SpO2 Anomaly	Temp Anomaly
08/01/2023 - 08/07/2023	1	1	0	0	0
08/08/2023 - 08/14/2023	0	0	0	0	0
08/15/2023 - 08/21/2023	0	0	0	0	0
08/22/2023 - 08/28/2023	0	0	0	0	0



Assessment stats on the Reporting tab will provide statistics from the results of clinically validated assessments that the patients have self-administered. These include the BAM R and IOP, the PHQ-9 and 15, GAD-7, and URICA. You can view these filtered by time, by patient, or by assessment type. Performance from admission to discharge can also be reviewed by type of assessment.


pretaa

- Events
- Patients
- Reporting**
  - Events
  - Assessment Stats
- Settings
- Feedback
- Contact Pretaa
- Version 1.1.0
- Logout

## Reporting: Assessments

Filter By: 1 Month All patients selected

Patients Overview	81 PHQ - 9	54 BAM-R	81 GAD-7	81 URICA	0 PHQ 15	0 BAM IOP
-------------------	---------------	-------------	-------------	-------------	-------------	--------------



pretaa

- Events
- Patients
- Reporting**
  - Events
  - Assessment Stats
- Settings
- Feedback
- Contact Pretaa
- Version 1.1.0
- Logout



Performance from admission to discharge

PHQ-9	↓ 55%	Average Score
GAD-7	↓ 49%	Average Score
URICA	↓ 23%	Readiness Score
BAM-R	↓ 99%	Use Factors
	↓ 39%	Risk Factor


pretaa

- Events
- Patients
- Reporting**
  - Events
  - Assessment Stats
- Settings
- Feedback
- Contact Pretaa
- Version 1.1.0
- Logout

## Reporting: Assessments

Filter By: 1 Month All patients selected

Patients Overview	81 PHQ - 9	54 BAM-R	81 GAD-7	81 URICA	0 PHQ 15	0 BAM IOP
-------------------	---------------	-------------	-------------	-------------	-------------	--------------




# Elevate your One-on-One's

Harnessing Pretaa's capabilities, clinicians can concentrate more on personalized patient care by reducing administrative burden, allowing for more engaging and meaningful one-on-ones.


Through real-time biometric data analysis and assessment scoring, Pretaa provides profound insights into patients' health status, empowering clinicians to tailor their interactions to address individual needs and concerns more effectively.


Additionally, Pretaa's proactive alert system enhances one-on-one sessions by identifying potential health risks, ensuring discussions are always relevant, timely, and focused on enhancing patients' recovery journey.



## Streamlining Patient Care for You

Designed to streamline administrative tasks for clinicians, optimizing your focus on patient care. Through technology, Pretaa proactively alerts clinicians about patient risk factors, enabling swift intervention.





## Finally, You Can Be Proactive not Reactive

Embrace proactive intervention with Pretaa, a sophisticated tool designed specifically for SUD clinicians, delivering alerts about possible crises before they escalate.

Leveraging biometric data, Pretaa accurately tracks parameters like Heart Rate, HR Variability, Skin Temp, SPO2, Steps, and Sleep, instantly notifying you about anomalies from the baseline 'normal'.

Further fortified with the latest in automated scoring technology for clinically validated assessments, Pretaa offers near-real-time alerts for significant risk factors, including suicidality indicators.

With Pretaa, you can shift from being reactive to proactive, enabling you to foresee complications, intervene promptly, and significantly improve patient outcomes.

## Self-Reporting Assessments Are Actually Self Reported

Pretaa revolutionizes the way SUD clinicians approach assessments by providing a platform for patients to self-administer clinically validated assessments such as BAM, PHQ-9 and URICA.

By leveraging Pretaa, clinicians can now bypass time-consuming EMR data entry during crucial treatment sessions and redirect their focus entirely on patient care and treatment strategies.

This modern, digital approach not only boosts efficiency but also empowers patients, allowing them to actively engage in their diagnostic process, paving the way for more personalized, effective treatment plans.





# Spend Less Time Diagnosing and More Time Treating.

We understand that as clinicians, your passion lies in helping your patients on their journey to recovery, not being buried in piles of paperwork.

From detailed protocols and assessments to the complexity of insurance forms, we know how these administrative tasks can distract you from what truly matters - the direct care of your patients.

We recognize that your higher calling is to heal and assist those battling addiction, and we're here to ensure that paperwork does not get in the way.



## Manage Your Caseload Smarter Not Harder

Discover the smarter way to manage your caseload with Pretaa's Clinicians Page, where Significant Patient Events and Detailed Reports offer critical insights into your patients' biometric data and clinically validated assessments.

Our platform empowers you to quickly identify trends within your caseload, ensuring you're always ahead in addressing patient needs and can swiftly act on high-risk situations.

By streamlining time management and enhancing focus, Pretaa transforms your valuable clinical hours into high-quality, effective patient care.

